



How to Install Version 6 in Windows Vista



It seems like everyone in the world is scrambling to make their programs work with the new Vista operating system from Microsoft. It looks like we are among the lucky ones. There are only a few things you need to do before installing Version 6 on a computer with the Microsoft Vista operating system.

Here is a summary of the things that need to be done.

- 1) Download and install the **.NET Framework 1.1 Service Pack 1** update to Vista (see link below)
- 2) Download and install **Microsoft SQL Server 2005 Express Edition Service Pack 2**
- 3) Temporarily turn off **User Account Control** in Vista.
- 4) Temporarily turn off all **Firewalls** and **Virus protection**.

Step 1: Download and install the **Microsoft .NET Framework 1.1 Service Pack 1** update which can be found by going to Microsoft's website, or for your convenience, we have provided a link to the file on our website.

1. Go to www.labsysgrp.com and point your mouse at **Support** and then click on **Support Documentation**.
2. Click on the link for the **.NET Framework 1.1 Service Pack 1**.
3. Install or **Run** the update. *You may need to restart your computer before the change will take place.*

Step 2: This update should already be on your computer if you have been keeping your Vista computer up to date with all of the latest Microsoft service patches. If you are not sure, then go to the Microsoft Vista Update site and let it update your computer with all of the latest patches and security packs. To check for updates:

1. Clicking on the **Start** button
2. Click on **Control Panel**.
3. In the Security area, click on **Check for updates**.

Step 3: Until Microsoft provides us with a better solution, you will need to turn off the Vista **User Account Control** before you will be able to install and use Version 6. To accomplish this, follow these steps:

- 1) Click on the **Start** button
- 2) Click on **Control Panel**
- 3) Click on **Security**
- 4) Click on **User Accounts** on far left side of the screen
- 5) Click on **User Accounts** on the right side of the screen
- 6) Click on the option to "**Turn User Accounts on or off**"
- 7) Remove check from box labeled "**Use User Account Control (UAC)**" to turn it off
- 8) Click **OK** to save and close the window.

Step 4: Temporarily turn off any and all Firewalls and virus protection software. These items can be turned back on after Lab Manager has been installed and is working properly. To turn off the Windows Firewall:

1. Click on the **Start** button
2. Click on **Control Panel**.
3. In the Security area, click on **Check the computers security status**.
4. If the Firewall is on, then click on **Windows Firewall** (on the left side of the screen).
5. Click on **Turn Windows Firewall on or off** (on the left side of the screen).
6. Select the **Off** option and click **OK**

Version 6 should now install and work properly in Windows Vista. If you continue to have problems, call our support department at 1-602-264-5913 and a support technician will help you.

